


KEOFITT REFLEX™ SAMPLING VALVE

USER MANUAL



INTRODUCTION:

MANUFACTURER: Keofitt A/S,
Lumbyvej 11H
5000 Odense C
Denmark

TYPE: KEOFITT REFLEX SAMPLING VALVE

PATENTS: EP1690038

YEAR OF INTRODUCTION: 2007

YEAR OF REVISED DESIGN: -

CONTENTS

Presentation.....	7
Valve function.....	7
Everyday use of the valve	8
Sterilisation	8
Sampling	8
Technical Data	9
Keofitt Reflex™ Pneumatic Valves	10
Parts and accessories for Keofitt Reflex™	12
Mounting instructions	13
Location:.....	13
Before welding:	13
Welding instructions	13
Maintenance	14
Spare parts list	14
Disassembly and assembly of valve body and head.....	14
Membranes	15

PRESENTATION

The Keofitt Reflex™ sampling valve is a unique 2-in-1 valve which can be readily sterilised and which meets both hygienic and process design requirements. This means that an effective cleaning and sterilisation of the sampling valve can be carried out between random samples independently of the course of the production process without compromising the same. The valve is 3-A and EHEDG Type EI authorised. 3-A Sanitary Standard is an American standard which is normative for a component's ease of cleaning and sterilisation. The standard ensures optimum conditions for food products which may come into contact with the component in question. The EHEDG Type EI certification is a European standard and it includes additional tests of bacterial increase on components that are in direct contact with the sample after the CIP process. The valve is used in a wide range of processing industries, such as breweries, dairies, and the pharmaceutical and biotechnological industries.

VALVE FUNCTION

The valve is designed to regularly take representative random samples in the production process. The valve is therefore designed such that effective cleaning, sterilisation and sampling can be carried out regularly without interrupting the production process.

Sterilisation is carried out by supplying steam through the steam inlet of the valve. It is the perfect, hygienic design and surface finish of the inner part of the valve which enables absolute sterilisation in a closed state. According to an EHEDG-based test carried out by the Biotechnological Institute in Denmark, the valve is sterilised after just 1 minute's supply of steam at a pressure of 2 bar(g) (121 °C).

Following sterilisation the valve is opened and the liquid will run out of the lower hose piece.

Note! The membrane functions both as a dynamic sample seal in the valve seat and as a hygienic, static steam seal against the valve body.



UNIQUE!

- The Keofitt Reflex valve controls both sample and steam supply. Also it eliminates the risk and waste involved with product going into the steam connection of other designs.



WARNING!

- During sterilisation with steam the valve will become hot, and care should thus be taken when handling the valve.
- The valve is designed for use in working conditions of up to 6 bar(g) pressure and temperatures of up to 121 °C. It is therefore important to be aware that the rubber plug (designed for max. 3 bar(g)) or the steel plug (designed for max. 10 bar(g)) can be forced out at high speed if not fitted correctly.
- Always remember to use safety goggles when taking samples because of the risk to the eyes.



IMPORTANT!

- The valve cannot be used for vacuum since the membrane will be sucked hard into the seat.
- The membranes are available in three different qualities: silicone, EPDM and PTFE.
- The silicone membrane has the advantage that it in general can stand higher temperatures, but it cannot tolerate moisture condensation resulting from steam sterilisation.
- The EPDM membrane is better able to cope with condensation in the steam, and at the same time can be used with a majority of CIP fluids.
- The PTFE membrane resists most CIP fluids and very high steam temperatures.

EVERYDAY USE OF THE VALVE

Sterilisation

Sterilisation takes place with valve in steam mode.

1. Open the steam supply and let it flow through the valve for sterilisation. 1 min. at 121 °C (2 bar(g)).
2. Close the steam supply.



WARNING!

- During sterilisation with steam the valve will become hot, and care should thus be taken when handling the valve.
- The valve is designed for use in working conditions of up to 6 bar(g) and temperatures of up to 121 °C. Always remember to use safety goggles when taking samples because of the risk to the eyes.
- Use saturated steam without condensation at max. 2 bar(g). At higher pressures the membrane can be damaged/split. In most cases the coaxial design ensures absolute cleanliness without the use of CIP or similar. If CIP is used, please refer to enclosed data sheet. If in doubt, contact Keofitt.




IMPORTANT!


- To reach 121°C a pressure of 2 bar(g) is needed. This can only be reached by use of a pressure release valve, or other counter pressure device mounted on the outlet connection.

Sampling

1. Open the valve and take the sample.
2. Shut the valve after the sample has been taken.
3. Clean the valve with steam and/or hot water, cf. 'sterilisation', points 1-4.



Sterilisation



Sampling

TECHNICAL DATA

Material

Valve body:	AISI 316L (1.4404)
Valve head:	AISI 316L (1.4404)
Membrane:	Silicone (grey) EPDM (black) PTFE (white)

Certificate

Valve body:	3.1
Membrane silicone	acc. to FDA & BGA
Membrane EPDM	acc. to FDA & BGA
Membrane PTFE	acc. to FDA & BGA



* A 6-digit code is marked on the valve body. This code refers to a 3.1 certificate which accompanies every consignment of valve bodies, certificates are available on www.keofitt.dk

Pressure - max.

Working pressure:	6 bar(g) / 87 psi(g)
Rubber plug	3 bar(g) / 44 psi(g)
Steel plug	10 bar(g) / 145 psi(g)

Temperature - max.

Sterilisation temp.: 121 °C / 249,80 °F **

** It is important that the steam is saturated, but dry, as condensation can damage the membrane. (Dry steam at max. 2 bar(g)).

Surface finish

Internal:	Ra ≤0,5 µm/20 µinch
External:	Electropolished


Every valve is Ra measured.

A serial number identifies each valve body. A surface finish certificate copy is available.

Viscosity:

Viscosity range: 0-1000cP, with particles up to 3mm in diameter.

Water at 20°C/68°F



$$K_v = Q \sqrt{\frac{p}{1000 * \Delta p}}$$

Kv: Valve capacity [m³/h]













Cv: Valve capacity [USgal/min]

Q: • Flow through valve seat [m³/h]


p: • Viscosity of fluid [kg/m³]

p: • Pressure drop across valve [bar]

KEOFITT REFLEX™ PNEUMATIC VALVES

		Hose piece				Mini tri-clamp				Welding liner			
Tank welding (type T) - process connection Silicone EPDM PTFE Key measurements <i>(Please refer to bottom page 11!)</i>													
	880001				880006				880007				
	880001E				880006E				880007E				
	885501				885501				885507				
	A	B	C	D	A	B	C	D	A	B	C	D	
14 mm	7 mm	38 mm	8 mm	25 mm	7 mm	38 mm	8 mm	12,7 mm	9,4 mm	38 mm	8 mm		
Pipe welding (type P) - process connection Silicone EPDM PTFE Key measurements <i>(Please refer to bottom page 11!)</i>													
	880011				880010				800017				
	880011E				880010E				880017E				
	885511				885510				885517				
	A	B	C	D	A	B	C	D	A	B	C	D	
14 mm	7 mm	38 mm	8 mm	25 mm	7 mm	38 mm	8 mm	12,7 mm	9,4 mm	38 mm	8 mm		
Clamp - connection 1" - process connection Silicone EPDM PTFE Key measurements <i>(Please refer to bottom page 11!)</i>													
	880021				880022				880027				
	880021E				880022E				880027E				
	885521				885522				885527				
	A	B	C	D	A	B	C	D	A	B	C	D	
14 mm	7 mm	50,4 mm	8 mm	25 mm	7 mm	50,4 mm	8 mm	12,7 mm	9,4 mm	50,4 mm	8 mm		
Clamp - connection 2" - process connection Silicone EPDM PTFE Key measurements <i>(Please refer to bottom page 11!)</i>													
	880003				880005				On request				
	880003E				880005E				On request				
	885503				885505				On request				
	A	B	C	D	A	B	C	D	A	B	C	D	
14 mm	7 mm	64 mm	8 mm	25 mm	7 mm	64 mm	8 mm	12,7 mm	9,4 mm	64 mm	8 mm		

	Hose piece	Mini tri-clamp	Welding liner																	
Varivent® connection Ø50 - process connection Silicone EPDM PTFE Key measurements <i>(Please refer to bottom page 11!)</i>																				
	880008																			
	880008E																			
	885508																			
	<table border="1"> <thead> <tr> <th>A</th> <th>B</th> <th>C</th> <th>D</th> </tr> </thead> <tbody> <tr> <td>14 mm</td> <td>7 mm</td> <td>50 mm</td> <td>8 mm</td> </tr> </tbody> </table>	A	B	C	D	14 mm	7 mm	50 mm	8 mm											
A	B	C	D																	
14 mm	7 mm	50 mm	8 mm																	
Varivent® connection Ø68 - process connection Silicone EPDM PTFE Key measurements <i>(Please refer to bottom page 11!)</i>																				
	880009																			
	880009E																			
	885509																			
	<table border="1"> <thead> <tr> <th>A</th> <th>B</th> <th>C</th> <th>D</th> </tr> </thead> <tbody> <tr> <td>14 mm</td> <td>7 mm</td> <td>68 mm</td> <td>8 mm</td> </tr> </tbody> </table>	A	B	C	D	14 mm	7 mm	68 mm	8 mm											
A	B	C	D																	
14 mm	7 mm	68 mm	8 mm																	
Ingold® connection Silicone EPDM PTFE Key measurements <i>(Please refer to bottom page 11!)</i>																				
		885532	885537																	
		<table border="1"> <thead> <tr> <th>A</th> <th>B</th> <th>C</th> <th>D</th> <th>A</th> <th>B</th> <th>C</th> <th>D</th> </tr> </thead> <tbody> <tr> <td>25 mm</td> <td>7 mm</td> <td>25 mm</td> <td>8 mm</td> <td>12,7 mm</td> <td>7 mm</td> <td>25 mm</td> <td>7 mm</td> </tr> </tbody> </table>	A	B	C	D	A	B	C	D	25 mm	7 mm	25 mm	8 mm	12,7 mm	7 mm	25 mm	7 mm		
	A	B	C	D	A	B	C	D												
25 mm	7 mm	25 mm	8 mm	12,7 mm	7 mm	25 mm	7 mm													




PARTS AND ACCESSORIES FOR KEOFITT REFLEX™

	Item	Material	Ident no.	Remark
	Membrane EPDM black	EPDM	880052	
	Membrane silicone grey	Silicone	880051	
	Membrane teflon	Teflon	880055	
	Rubber cap with knob		600062	
	Adjustable pressure relief valve (0-3 Bar(g))		850059	
	Quick coupling steel pipe		800070	
	Quick coupling for PTFE		800071	
	Quick coupling with 3/8" Hose Barb		800082	
	Sampling coil with st.st. cap		800058	
	Quick coupling with M4 hose piece		800086	
	O-ring	Silicone	600825	
	O-ring	Silicone	880830	

MOUNTING INSTRUCTIONS

Location:

The valve should always be located with its centre line in a horizontal or lower position. The valve will then be selfdraining.



Before welding:

Remember to disassemble the valve body and head. The valve body and head must be separated during welding. Rubber plugs, chain and membrane must be removed from the valve body, as otherwise heat from the welding process will damage them.

WELDING INSTRUCTIONS

Valves for welding are available in two types: T (tank) and P (pipe).

1. For type T (tank) it is necessary to drill a hole $\varnothing 28$ mm into the tank wall, and then fit the valve into this hole flush with the inside of the tank. Welding should be carried out as a penetration welding.
Material thickness less than 4 mm: Weld from inside. Material thickness greater than 4 mm: Weld from both outside and inside.
Since type T has a solid end piece, the valve will not be damaged by penetration welding. However, the use of purge gas in the form of either Argon or Formier gas is recommended in order to give the best result.
2. For type P (pipe) penetration welding must be carried out from outside. The valve is machined with a recess-like shoulder on the outside of the end piece which gives approximately the same material thickness (1.5mm material thickness) as in the pipe wall.




IMPORTANT!

- When grinding/polishing the internal weld, the valve seat must not be touched.

The welding result will be best if the following method is used:

A collar is made on the pipe section so that the valve has a flat contact face. This flaring must look like a T-piece, as shown in the example below.




- The pipe section and the valve's hose pieces are sealed with sponge rubber or similar.
- Purge gas such as Argon or Formier gas is fed through the valve body into the pipe section and the system is now filled with 6 times the estimated volume of the pipe section. All O₂ is thus expelled from the system and welding can commence.
- Welding can take place with the purge gas continually flowing in the system.
- The gas remains in the system until the item is lukewarm, after which the set-up can be dismantled.

MAINTENANCE

The rubber membrane should be replaced every two months. PTFE membranes should be replaced every 12 months. In the event of intensive sterilisation and cleaning it may be necessary to replace it more frequently. The rubber plug must be replaced at least once every six months. In each individual case a standard operating procedure including maintenance intervals should be endorsed based on experience. For disassembly of valve body and valve head, see instructions.


Spare parts



Pos. Item

- 1. O-ring
- 2. Membrane Silicone (grey)
Membrane EPDM (black)
Membrane PTFE (White)

Disassembly and assembly of valve body and head



ORDER OF OPERATION:

Remember! When replacing the membrane, set the valve head in the open position before disassembling from the body. Omitting to do so may result in twisting and cutting of the membrane.

- 1. Set the valve head at the open position.
- 2. Remove the valve head pos. 1. This is carried out by loosening the clamp and then pulling the valve head off.
- 3. Refit the valve head (in the open position) once the necessary parts have been replaced.

SILICONE MEMBRANE - ITEM NO. 880051




TECHNICAL SPECIFICATION

- | | |
|--|-----------------------|
| • Type: | Silicone |
| • Colour: | Grey |
| • Resistance to chemicals acids/bases) | Suitable |
| • Food safe | Yes (FDA*) |
| • Temp. max. | 130-250°C / 266-482°F |
| • Steam pressure max. | 2 bar/29 psi |
| • Process pressure | 1-6 bar/14-87 psi |

Keofitt recommends to change the silicone membrane 4-6 times a year or as needed. The recommendation is based on 1-5 samples a day, but should reflect individual tear and wear from individual cleaning and sterilisation procedures.

*FDA approved compound according to Code of Federal Regulations Title 21 - § 177.1550



EPDM MEMBRANE - ITEM NO. 880052




TECHNICAL SPECIFICATION

- | | |
|--|-----------------------|
| • Type: | EPDM |
| • Colour: | Black |
| • Resistance to chemicals acids/bases) | Very good |
| • Food safe | Yes (FDA*) |
| • Temp. max. | 130-150°C / 266-302°F |
| • Steam pressure max. | 2 bar/29 psi |
| • Process pressure | 1-6 bar/14-87 psi |

Keofitt recommends to change the EPDM membrane 4-6 times a year or as needed. The recommendation is based on 1-5 samples a day, but should reflect individual tear and wear from individual cleaning and sterilisation procedures.


*FDA approved compound according to Code of Federal Regulations Title 21 - § 177.1550



PTFE MEMBRANE - ITEM NO. 880055

TECHNICAL SPECIFICATION

- | | |
|--|-----------------------|
| • Type: | PTFE |
| • Colour: | White |
| • Resistance to chemicals acids/bases) | Excellent** |
| • Food safe | Yes (FDA*) |
| • Temp. max. | 130-250°C / 266-482°F |
| • Steam pressure max. | 2 bar/29 psi |
| • Process pressure | 1-6 bar/14-87 psi |




Keofitt recommends to change the PTFE membrane once a year or as needed. The recommendation is based on 1-5 samples a day, but should reflect individual tear and wear from individual cleaning and sterilisation procedures.

*FDA approved compound according to Code of Federal Regulations Title 21 - § 177.1550

** Is not attacked by common chemicals, with the exception of strongly oxidising acids.

INSTRUCTIONS ON REPLACING PTFE MEMBRANE

1. Open valve.
2. Remove the valve head from the valve body.
3. Push the membrane upwards until the tool for membrane fits under it.
4. Insert tool for membrane, between the membrane and the bushing.
5. Close valve head.
6. Now the membrane is loosened from the valve head and can be replaced.



To attach new membrane to valve head.

7. Set the valve head to closed position.
8. Place the new membrane on valve head.
9. Mount the membrane bushing with the new Teflon membrane by pressing the membrane with your hand until it clicks.
10. Set the valve head in open position.
11. Insert the valve head into the valve body.
12. Close valve head.

IMPORTANT

- Once the membrane has been removed from the valve head the click system in the membrane might be damaged. Therefore the membrane might be unsafe for further use and it is recommended not to use the membrane again.
- Do not use hammer or other tool that might scratch the surface of the membrane.



Tool for membrane item no. 600255

Keofitt reserves the right to change technical data without notice!
For complete set of updated data sheets and manuals for Keofitt products please refer to our web page www.keofitt.dk

